



# NanoPhotometer<sup>®</sup> NP80

*distributed by*



## All-in-one UV-Vis Spectroscopy for Quality Control in the Wine and Spirits Industry

*Confirm product consistency and identify deviations  
during various stages of production*

- ▶ **Higher Precision:** No dilutions needed / No calibration / Only small amount of liquid used
- ▶ **Save Time and Money:** One time investment / No consumables / Maintenance free



# NanoPhotometer<sup>®</sup> NP80

## All-in-one UV-Vis Spectroscopy for Quality Control in the Wine and Spirits Industry

- No dilutions required
- Rapid and complete sample analysis in as little as 2.5 seconds
- NanoVolume (5µl) and cuvette functionality for sample analysis
- Scan range: 200 – 900 nm
- Zero cost of ownership: No recalibration or regular maintenance required – ever
- World's smallest footprint in its class: 20 x 20 x 12 cm
- Stand-alone operation with built-in 7-inch glove compatible touch screen
- Mobile use (battery powered)
- Endless Connectivity (WiFi, HotSpot, LAN)
- Built-in vortex to achieve sample uniformity



## UV/Vis Spectrophotometry

The NanoPhotometer NP80 has a wavelength range of 200 to 900 nm, and is suited to perform standard UV/Vis photometric (quantitation of anthocyanins, phenols, flavonoids, Folin Ciocalteu etc.), colorimetric (color parameters such as color pigments, color density, intensity/hue) and enzymatic (Lactic Acid, Ethanol, Glucose, Fructose, Sucrose, etc.) analytical procedures common in wine, spirit, and juice analysis.

The available methods include wavelength (fixed wavelengths), wavescan, kinetics and standard curve applications and all methods can be customized by the user to be compatible with most ready-to-use commercial reagent kits or for method development. Patented Sample Compression Technology™ provides unmatched precision and accuracy. A capillary film is formed between two quartz surfaces with each sample, eliminating the need for surface tension and avoiding evaporation.



# ImplenQ Software

ImplenQ is a cutting-edge software package specifically designed for quality control (QC) laboratories in a wide range of industries. Its primary function is to ensure the consistency of products by leveraging UV/Vis spectrophotometry for verifying batch-to-batch quality and identifying deviations in samples through spectral fingerprinting. A sample is scanned over the entire wavelength range from 200–900 nm to capture its unique spectral signature within 2.5 seconds. ImplnQ is designed for NanoVolume measurements, requiring a small droplet of only 5µl, without the need for further dilution steps.



## Recognition Algorithm



The ImplnQ software employs algorithms to compare the generated UV/Vis spectral fingerprints to the database of reference spectra. The algorithm is able to identify microscopic deviations in the spectral fingerprints, thereby providing the user invaluable insights into their product composition.



## Data Management

ImplenQ stores and organizes reference spectra in a user-friendly database. Each reference spectrum serves as a digital fingerprint for a specific product sample. This database serves as a repository of trusted information for comparative analyses.

## Example of Product Comparisons



Comparison of spectral fingerprint scans performed on both authentic and counterfeit wine (A) and whiskey (B) samples. Samples were scanned over the entire wavelength range of 200–900 nm and the spectral fingerprints were analyzed by the ImplnQ software algorithm to identify the counterfeit products.

# Technical Specifications

NanoVolume Performance		Cuvette Performance	
Sample Volume	5 µl	Photometric Range	0 - 2.6 A
Photometric Range (10 mm equivalent)	0.02 - 330 A	Center Height (Z-Height)	8.5 mm
Path Length	0.67 & 0.07 mm	Cell Types	Outside dimension 12.5 x 12.5 mm
Dilution Factor	15 and 140	Heating	37 °C ± 0.5 °C
Vortex	2,800 rpm Tube size up to 2.0 ml	Processing Power & Compatibility	
Optical Specifications		Operating System	Linux based NPOS
Wavelength Scan Range	200 - 900 nm	Onboard Processor	Intel Celeron dual core 2.4 GHz
Measure Time For Full Scan Range	2.5 - 4.0 sec	Internal Data Storage	64 GB
Wavelength Reproducibility	± 0.2 nm	In & Output Ports	2x USB A, USB B, HDMI, Ethernet, WiFi
Wavelength Accuracy	± 0.75 nm	Software Compatibility	Windows 8, 10 (32 & 64 bit) and 11 OS X (Intel x86 and Apple M1)
Bandwidth	< 1.5 nm	General Specifications	
Absorbance Reproducibility	(Cuvette): < 0.002 A @ 0 - 0.3 A @ 280 nm CV < 1% @ 0.3 - 2.0 A @ 280 nm  (Lid 15): < 0.002 A @ 0 - 0.3 A @ 280 nm CV < 1% @ 0.3 - 1.7 A @ 280 nm	Main Body Size	200 x 200 x 120 mm
Absorbance Accuracy	< 1.75% @ 0.7 A @ 280 nm of the reading	Weight	3.8 - 5.2 kg depending on configuration
Stray Light	< 0.5% @ 240 nm using NaI	Operating Voltage	90 - 250 V ± 10%, 50/60 Hz, 90 W, 18/19 VDC
Optical Arrangement	1x 4096 CMOS Array	Display	1024 x 600 pixels; glove compatible touchscreen
Lamp   Lifetime	Xenon flash lamp 10 <sup>9</sup> flashes, up to 10 years	Built-in Battery Pack: Optional rechargeable lithium ion battery	95 Wh, 6.6 Ah, 8 h Min. charging cycles: 800
		Certification	CE, IEC 61010-1:2012 and EN 61326-1:2013
		Battery Certification	IEC 62133 and UN38.3 transport test
		Security	Slot for Kensington lock



## NanoPhotometer<sup>®</sup> NP80

*distributed by*



Contact us today!

info@bevzero.com | (707) 577-7500 | BEVZERO.COM